

产品名称: Cyclin A2 Rabbit Monoclonal Antibody
产品货号: AMRe86267



产品概述 (Summary)

产品名称 (Production Name)	Cyclin A2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cyclin A2
别名 (Alternative Names)	CCN1; CCNA
基因 ID (Gene ID)	890
蛋白 ID (SwissProt ID)	P20248.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:2000,ICC/IF 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:49 kDa; Observed MW:49 kDa

研究背景 (Background)


产品名称: Cyclin A2 Rabbit Monoclonal Antibody
产品货号: AMRe86267




The protein encoded by this gene belongs to the highly conserved cyclin family, whose members function as regulators of the cell cycle. This protein binds and activates cyclin-dependent kinase 2 and thus promotes transition through G1/S and G2/M. [provided by RefSeq, Aug 2016]

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe86267 at 1:1000.



Western blot analysis of extracts from T24 cells using AMRe86267 at 1:1000.



Immunofluorescence analysis of HeLa cells labelling Cyclin A2 with AMRe86267.

产品名称: Cyclin A2 Rabbit Monoclonal Antibody
产品货号: AMRe86267



注意事项 (Note)

For research use only .