

产品名称: SAM68 Rabbit Monoclonal Antibody
产品货号: AMRe86275



产品概述 (Summary)

产品名称 (Production Name)	SAM68 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SAM68
别名 (Alternative Names)	p62; p68; Sam68
基因 ID (Gene ID)	10657
蛋白 ID (SwissProt ID)	Q07666.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:5000-1:50000,IHC 1:100-1:200,ICC/IF 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:49 kDa; Observed MW:68 kDa

研究背景 (Background)


产品名称: SAM68 Rabbit Monoclonal Antibody
产品货号: AMRe86275



This gene encodes a member of the K homology domain-containing, RNA-binding, signal transduction-associated protein family. The encoded protein appears to have many functions and may be involved in a variety of cellular processes, including alternative splicing, cell cycle regulation, RNA 3'-end formation, tumorigenesis, and regulation of human immunodeficiency virus gene expression. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2012]

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of HeLa cells labelling SAM68 with AMRe86275.

注意事项 (Note)

For research use only .