

产品名称: HAS3 Rabbit Monoclonal Antibody
产品货号: AMRe86285



产品概述 (Summary)

产品名称 (Production Name)	HAS3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HAS3
别名 (Alternative Names)	HA synthase 3
基因 ID (Gene ID)	3038
蛋白 ID (SwissProt ID)	O00219.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW:63 kDa; Observed MW:63 kDa

研究背景 (Background)


产品名称: HAS3 Rabbit Monoclonal Antibody
产品货号: AMRe86285



The protein encoded by this gene is involved in the synthesis of the unbranched glycosaminoglycan hyaluronan, or hyaluronic acid, which is a major constituent of the extracellular matrix. This gene is a member of the NODC/HAS gene family. Compared to the proteins encoded by other members of this gene family, this protein appears to be more of a regulator of hyaluronan synthesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from PC-3 cells using AMRe86285 at 1:1000.

注意事项 (Note)

For research use only .