

产品名称: PODXL Rabbit Monoclonal Antibody
产品货号: AMRe86294



产品概述 (Summary)

产品名称 (Production Name)	PODXL Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PODXL
别名 (Alternative Names)	PC; PCLP; Gp200; PCLP-1
基因 ID (Gene ID)	5420
蛋白 ID (SwissProt ID)	O00592.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:1000,ICC/IF 1:100-1:200,FC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:59 kDa; Observed MW:160 kDa

研究背景 (Background)


产品名称: PODXL Rabbit Monoclonal Antibody
产品货号: AMRe86294



This gene encodes a member of the sialomucin protein family. The encoded protein was originally identified as an important component of glomerular podocytes. Podocytes are highly differentiated epithelial cells with interdigitating foot processes covering the outer aspect of the glomerular basement membrane. Other biological activities of the encoded protein include: binding in a membrane protein complex with Na⁺/H⁺ exchanger regulatory factor to intracellular cytoskeletal elements, playing a role in hematopoietic cell differentiation, and being expressed in vascular endothelium cells and binding to L-selectin. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of HeLa cells labelling PODXL with AMRe86294.

注意事项 (Note)

For research use only .