

产品名称: MAGEA3 Rabbit Monoclonal Antibody
产品货号: AMRe86300



产品概述 (Summary)

产品名称 (Production Name)	MAGEA3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAGEA3
别名 (Alternative Names)	HIP8; HYPD; CT1.3; MAGE3; MAGEA6
基因 ID (Gene ID)	4102
蛋白 ID (SwissProt ID)	P43357.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:500,ICC/IF 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:35 kDa; Observed MW:45 kDa

研究背景 (Background)


产品名称: MAGEA3 Rabbit Monoclonal Antibody
产品货号: AMRe86300




This gene is a member of the MAGEA gene family. The members of this family encode proteins with 50 to 80% sequence identity to each other. The promoters and first exons of the MAGEA genes show considerable variability, suggesting that the existence of this gene family enables the same function to be expressed under different transcriptional controls. The MAGEA genes are clustered at chromosomal location Xq28. They have been implicated in some hereditary disorders, such as dyskeratosis congenita. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from A375 cells using AMRe86300 at 1:1000.



Immunofluorescence analysis of HeLa cells labelling MAGEA3 with AMRe86300.

注意事项 (Note)

For research use only .