

产品名称: DNA PKcs Rabbit Monoclonal Antibody
产品货号: AMRe86303



产品概述 (Summary)

产品名称 (Production Name)	DNA PKcs Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DNA PKcs
别名 (Alternative Names)	HYRC; p350; DNAPK; DNPk1; HYRC1; IMD26; XRCC7; DNA-PKcs
基因 ID (Gene ID)	5591
蛋白 ID (SwissProt ID)	P78527.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100,ICC/IF 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:469 kDa; Observed MW:469 kDa

研究背景 (Background)


产品名称: DNA PKcs Rabbit Monoclonal Antibody
产品货号: AMRe86303



This gene encodes the catalytic subunit of the DNA-dependent protein kinase (DNA-PK). It functions with the Ku70/Ku80 heterodimer protein in DNA double strand break repair and recombination. The protein encoded is a member of the PI3/PI4-kinase family.[provided by RefSeq, Jul 2010]

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of HeLa cells labelling DNA PKcs with AMRe86303.

注意事项 (Note)

For research use only .