

产品名称: Lin28A Rabbit Monoclonal Antibody
产品货号: AMRe86310



产品概述 (Summary)

产品名称 (Production Name)	Lin28A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Lin28A
别名 (Alternative Names)	CSDD1; LIN28; LIN-28; ZCCHC1; lin-28A
基因 ID (Gene ID)	79727
蛋白 ID (SwissProt ID)	Q9H9Z2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000, ICC/IF 1:100-1:500, FC 1:10-1:100, IP 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:23 kDa; Observed MW:26 kDa

研究背景 (Background)


产品名称: Lin28A Rabbit Monoclonal Antibody
产品货号: AMRe86310



This gene encodes a LIN-28 family RNA-binding protein that acts as a posttranscriptional regulator of genes involved in developmental timing and self-renewal in embryonic stem cells. The encoded protein functions through direct interaction with target mRNAs and by disrupting the maturation of certain miRNAs involved in embryonic development. This protein prevents the terminal processing of the LET7 family of microRNAs which are major regulators of cellular growth and differentiation. Aberrant expression of this gene is associated with cancer progression in multiple tissues. [provided by RefSeq, Sep 2015]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from NCCIT cells using AMRe86310 at 1:1000.

注意事项 (Note)

For research use only .