

产品名称: GRP94 Rabbit Monoclonal Antibody
产品货号: AMRe86314



产品概述 (Summary)

产品名称 (Production Name)	GRP94 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GRP94
别名 (Alternative Names)	ECGP; GP96; TRA1; GRP94; HEL35; HEL-S-125m
基因 ID (Gene ID)	7184
蛋白 ID (SwissProt ID)	P14625.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:93 kDa; Observed MW:94 kDa

研究背景 (Background)


产品名称: GRP94 Rabbit Monoclonal Antibody
产品货号: AMRe86314



This gene encodes a member of a family of adenosine triphosphate(ATP)-metabolizing molecular chaperones with roles in stabilizing and folding other proteins. The encoded protein is localized to melanosomes and the endoplasmic reticulum. Expression of this protein is associated with a variety of pathogenic states, including tumor formation. There is a microRNA gene located within the 5' exon of this gene. There are pseudogenes for this gene on chromosomes 1 and 15. [provided by RefSeq, Aug 2012]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from T47D cells using AMRe86314 at 1:1000.

注意事项 (Note)

For research use only .