

产品名称: CD79b Rabbit Monoclonal Antibody  
产品货号: AMRe86320



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | CD79b Rabbit Monoclonal Antibody       |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB, ICC/IF, FC                         |
| 种属反应性 (Reactivity)     | Human                                  |

## 产品性能 (Performance)

|                     |                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                                                                                   |
| 形式 (Form)           | Liquid                                                                                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                                                     |
| 储存溶液 (Buffer)       | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt. |
| 纯化方式 (Purification) | Affinity Purification                                                                                                                                        |

## 免疫原信息 (Immunogen)

|                        |                |
|------------------------|----------------|
| 基因名 (Gene Name)        | CD79b          |
| 别名 (Alternative Names) | B29; IGB; AGM6 |
| 基因 ID (Gene ID)        | 974            |
| 蛋白 ID (SwissProt ID)   | P40259.        |

## 产品应用 (Application)

|                          |                                                     |
|--------------------------|-----------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:5000, ICC/IF 1:200-1:500, FC 1:10-1:100 |
| 蛋白分子量 (Molecular Weight) | Calculated MW:26 kDa; Observed MW:30-50 kDa         |

## 研究背景 (Background)


**产品名称: CD79b Rabbit Monoclonal Antibody**  
**产品货号: AMRe86320**



The B lymphocyte antigen receptor is a multimeric complex that includes the antigen-specific component, surface immunoglobulin (Ig). Surface Ig non-covalently associates with two other proteins, Ig-alpha and Ig-beta, which are necessary for expression and function of the B-cell antigen receptor. This gene encodes the Ig-beta protein of the B-cell antigen component. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of extracts from Daudi cells using AMRe86320 at 1:1000.

## 注意事项 (Note)

For research use only .