

产品名称: ASMT Rabbit Monoclonal Antibody
产品货号: AMRe86335



产品概述 (Summary)

产品名称 (Production Name)	ASMT Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ASMT
别名 (Alternative Names)	ASMTY; HIOMT; HIOMTY
基因 ID (Gene ID)	438
蛋白 ID (SwissProt ID)	P46597.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:39 kDa; Observed MW:39 kDa

研究背景 (Background)


产品名称: ASMT Rabbit Monoclonal Antibody
产品货号: AMRe86335



This gene belongs to the methyltransferase superfamily, and is located in the pseudoautosomal region (PAR) at the end of the short arms of the X and Y chromosomes. The encoded enzyme catalyzes the final reaction in the synthesis of melatonin, and is abundant in the pineal gland. Alternatively spliced transcript variants have been noted for this gene. [provided by RefSeq, Jan 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Y79 cells using AMRe86335 at 1:1000.

注意事项 (Note)

For research use only .