

产品名称: GRSF1 Rabbit Monoclonal Antibody
产品货号: AMRe86348



产品概述 (Summary)

产品名称 (Production Name)	GRSF1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GRSF1
别名 (Alternative Names)	G-rich sequence factor 1
基因 ID (Gene ID)	2926
蛋白 ID (SwissProt ID)	Q12849.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:100-1:200, FC 1:50-1:200, IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:53 kDa; Observed MW:53 kDa

研究背景 (Background)


产品名称: GRSF1 Rabbit Monoclonal Antibody
产品货号: AMRe86348



The protein encoded by this gene is a cellular protein that binds RNAs containing the G-rich element. The protein is localized in the cytoplasm, and has been shown to stimulate translation of viral mRNAs in vitro. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of HeLa cells labelling GRSF1 with AMRe86348.

注意事项 (Note)

For research use only .