

产品名称: GFP Rabbit Monoclonal Antibody
产品货号: AMRe86353



产品概述 (Summary)

产品名称 (Production Name)	GFP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Other

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GFP
别名 (Alternative Names)	GFP
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P42212.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:10000,IHC 1:100-1:200,ICC/IF 1:200-1:500,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:27 kDa; Observed MW:According to the fused protein

研究背景 (Background)


产品名称: GFP Rabbit Monoclonal Antibody
产品货号: AMRe86353



Energy-transfer acceptor. Its role is to transduce the blue chemiluminescence of the protein aequorin into green fluorescent light by energy transfer. Fluoresces in vivo upon receiving energy from the Ca^{2+} -activated photoprotein aequorin.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of GFP in CHO-K1, 944-GFP lysates using GFP Tag antibody.

注意事项 (Note)

For research use only .