

产品名称: Lamin B1 Rabbit Monoclonal Antibody
产品货号: AMRe86356



产品概述 (Summary)

产品名称 (Production Name)	Lamin B1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Lamin B1
别名 (Alternative Names)	LMN; ADLD; LMN2; LMNB
基因 ID (Gene ID)	4001
蛋白 ID (SwissProt ID)	P20700.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:20-1:500,ICC/IF 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:66 kDa; Observed MW:70 kDa

研究背景 (Background)


产品名称: Lamin B1 Rabbit Monoclonal Antibody
产品货号: AMRe86356




This gene encodes one of the two B-type lamin proteins and is a component of the nuclear lamina. A duplication of this gene is associated with autosomal dominant adult-onset leukodystrophy (ADLD). Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2015]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of Lamin B1 in K562, Rat Brain, C6, 3T3, HeLa cell lysates using Lamin B1 antibody (1:1000 diluted).



Western blot analysis of extracts from HepG2 cells and Rat heart, Mouse kidney tissue using AMRe86356 at 1:1000.

注意事项 (Note)

For research use only .