

产品名称: S Tag Rabbit Monoclonal Antibody
产品货号: AMRe86358



产品概述 (Summary)

产品名称 (Production Name)	S Tag Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Transfected

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	S Tag
别名 (Alternative Names)	S Tag
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW;; Observed MW:According to the fused protein

研究背景 (Background)


产品名称: S Tag Rabbit Monoclonal Antibody
产品货号: AMRe86358



According to the fused protein

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of S tag in multi-tag antigen 120ng, multi-tag antigen 210ng cell lysates using S tag antibody.

注意事项 (Note)

For research use only .