

产品名称: Transferrin Rabbit Monoclonal Antibody
产品货号: AMRe86359



产品概述 (Summary)

产品名称 (Production Name)	Transferrin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Transferrin
别名 (Alternative Names)	TFQTL1; PRO1557; PRO2086; HEL-S-71p
基因 ID (Gene ID)	7018
蛋白 ID (SwissProt ID)	P02787.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:500,FC 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:77 kDa; Observed MW:77 kDa

研究背景 (Background)


产品名称: Transferrin Rabbit Monoclonal Antibody
产品货号: AMRe86359




This gene encodes a glycoprotein with an approximate molecular weight of 76.5 kDa. It is thought to have been created as a result of an ancient gene duplication event that led to generation of homologous C and N-terminal domains each of which binds one ion of ferric iron. The function of this protein is to transport iron from the intestine, reticuloendothelial system, and liver parenchymal cells to all proliferating cells in the body. This protein may also have a physiologic role as granulocyte/pollen-binding protein (GPBP) involved in the removal of certain organic matter and allergens from serum. [provided by RefSeq, Sep 2009]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of Transferrin in K562 cell lysates using Transferrin antibody(1:1000 diluted).



Western blot analysis of extracts from Human plasma tissue using AMRe86359 at 1:1000.

注意事项 (Note)

For research use only .