

产品名称: YY1 Rabbit Monoclonal Antibody
产品货号: AMRe86362



产品概述 (Summary)

产品名称 (Production Name)	YY1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	YY1
别名 (Alternative Names)	DELTA; NF-E1; UCRBP; GADEV5; INO80S; YIN-YANG-1
基因 ID (Gene ID)	7528
蛋白 ID (SwissProt ID)	P25490.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:200,ICC/IF 1:20-1:50,FC 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:45 kDa; Observed MW:65 kDa

研究背景 (Background)


产品名称: YY1 Rabbit Monoclonal Antibody
产品货号: AMRe86362




YY1 is a ubiquitously distributed transcription factor belonging to the GLI-Kruppel class of zinc finger proteins. The protein is involved in repressing and activating a diverse number of promoters. YY1 may direct histone deacetylases and histone acetyltransferases to a promoter in order to activate or repress the promoter, thus implicating histone modification in the function of YY1. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of YY1 in K562,Rat Brain,C6,3T3,HeLa cell lysates using YY1 antibody(1:1000 diluted).



Western blot analysis of extracts from T24 cells and Mouse heart tissue using AMRe86362 at 1:2000.

注意事项 (Note)

For research use only .