

产品名称: CDX2 Rabbit Monoclonal Antibody
产品货号: AMRe86371



产品概述 (Summary)

产品名称 (Production Name)	CDX2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDX2
别名 (Alternative Names)	CDX3; CDX-3; CDX2/AS
基因 ID (Gene ID)	1045, 12591
蛋白 ID (SwissProt ID)	Q99626, P43241.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:1000,ICC/IF 1:200-1:500
蛋白分子量 (Molecular Weight)	Calculated MW:34 kDa; Observed MW:38 kDa

研究背景 (Background)


产品名称: CDX2 Rabbit Monoclonal Antibody
产品货号: AMRe86371



This gene is a member of the caudal-related homeobox transcription factor gene family. The encoded protein is a major regulator of intestine-specific genes involved in cell growth and differentiation. This protein also plays a role in early embryonic development of the intestinal tract. Aberrant expression of this gene is associated with intestinal inflammation and tumorigenesis. [provided by RefSeq, Jan 2012]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Caco-2 cells using AMRe86371 at 1:1000.

注意事项 (Note)

For research use only .