

产品名称: SENP1 Rabbit Monoclonal Antibody
产品货号: AMRe86379



产品概述 (Summary)

产品名称 (Production Name)	SENP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SENP1
别名 (Alternative Names)	SuPr-2
基因 ID (Gene ID)	29843
蛋白 ID (SwissProt ID)	Q9P0U3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:10000,IHC 1:100-1:200,ICC/IF 1:100-1:500,FC 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:74 kDa; Observed MW:74 kDa

研究背景 (Background)


产品名称: SENP1 Rabbit Monoclonal Antibody
产品货号: AMRe86379




This gene encodes a cysteine protease that specifically targets members of the small ubiquitin-like modifier (SUMO) protein family. This protease regulates SUMO pathways by deconjugating sumoylated proteins. This protease also functions to process the precursor SUMO proteins into their mature form. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Jun 2012]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe86379 at 1:1000.



Immunofluorescence analysis of HeLa cells labelling SENP1 with AMRe86379.

注意事项 (Note)

For research use only .