

产品名称: Importin4 Rabbit Monoclonal Antibody
产品货号: AMRe86384



产品概述 (Summary)

产品名称 (Production Name)	Importin4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Importin4
别名 (Alternative Names)	Imp4
基因 ID (Gene ID)	79711
蛋白 ID (SwissProt ID)	Q8TEX9.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:50-1:100,ICC/IF 1:50-1:100,FC 1:10-1:100,IP 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:119 kDa; Observed MW:130 kDa

产品名称: Importin4 Rabbit Monoclonal Antibody
产品货号: AMRe86384




研究背景 (Background)

Predicted to enable nuclear import signal receptor activity and nuclear localization sequence binding activity. Involved in DNA replication-dependent chromatin assembly; DNA replication-independent chromatin assembly; and protein import into nucleus. Located in chromatin. Part of protein-containing complex. [provided by Alliance of Genome Resources, Apr 2022]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HepG2 cells using AMRe86384 at 1:1000.

注意事项 (Note)

For research use only .