

产品名称: MAD1 Rabbit Monoclonal Antibody
产品货号: AMRe86394



产品概述 (Summary)

产品名称 (Production Name)	MAD1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAD1
别名 (Alternative Names)	MAD1; PIG9; TP53I9; TXBP181
基因 ID (Gene ID)	8379
蛋白 ID (SwissProt ID)	Q9Y6D9.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:83 kDa; Observed MW:83 kDa

研究背景 (Background)


产品名称: MAD1 Rabbit Monoclonal Antibody
产品货号: AMRe86394



MAD1L1 is a component of the mitotic spindle-assembly checkpoint that prevents the onset of anaphase until all chromosome are properly aligned at the metaphase plate. MAD1L1 functions as a homodimer and interacts with MAD2L1. MAD1L1 may play a role in cell cycle control and tumor suppression. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2015]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HepG2, MCF-7, A431, HeLa cells using AMRe86394 at 1:1000.

注意事项 (Note)

For research use only .