

产品名称: Tartrate Resistant Acid Phosphatase Rabbit Monoclonal Antibody
产品货号: AMRe86395



产品概述 (Summary)

产品名称 (Production Name)	Tartrate Resistant Acid Phosphatase Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Tartrate Resistant Acid Phosphatase
别名 (Alternative Names)	HPAP; TRAP; TRACP5a; TRACP5b; TrATPase
基因 ID (Gene ID)	54
蛋白 ID (SwissProt ID)	P13686.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:200,ICC/IF 1:100-1:200,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:37 kDa; Observed MW:35 kDa

**产品名称: Tartrate Resistant Acid Phosphatase Rabbit
Monoclonal Antibody**
产品货号: AMRe86395




研究背景 (Background)

This gene encodes an iron containing glycoprotein which catalyzes the conversion of orthophosphoric monoester to alcohol and orthophosphate. It is the most basic of the acid phosphatases and is the only form not inhibited by L(+)-tartrate. [provided by RefSeq, Aug 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from A431 cells using AMRe86395 at 1:1000.

注意事项 (Note)

For research use only .