

产品名称: PCGF5 Rabbit Monoclonal Antibody
产品货号: AMRe86401



产品概述 (Summary)

产品名称 (Production Name)	PCGF5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PCGF5
别名 (Alternative Names)	AI324127; 0610009F02Rik; 1110054A01Rik; 5830406C17Rik; 5830443C21Rik; 9530023M17Rik
基因 ID (Gene ID)	76073
蛋白 ID (SwissProt ID)	Q3UK78.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:30 kDa; Observed MW:30 kDa

产品名称: PCGF5 Rabbit Monoclonal Antibody
产品货号: AMRe86401




研究背景 (Background)

Component of a Polycomb group (PcG) multiprotein PRC1-like complex, a complex class required to maintain the transcriptionally repressive state of many genes, including Hox genes, throughout development. PcG PRC1 complex acts via chromatin remodeling and modification of histones; it mediates monoubiquitination of histone H2A 'Lys-119', rendering chromatin heritably changed in its expressibility (By similarity).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from F9 cells using AMRe86401 at 1:1000.

注意事项 (Note)

For research use only .