

产品名称: DPP3 Rabbit Monoclonal Antibody
产品货号: AMRe86408



产品概述 (Summary)

产品名称 (Production Name)	DPP3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DPP3
别名 (Alternative Names)	DPPIII
基因 ID (Gene ID)	10072, 75221, 114591
蛋白 ID (SwissProt ID)	Q9NY33, Q99KK7, O55096.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW:83 kDa; Observed MW:83 kDa

研究背景 (Background)


产品名称: DPP3 Rabbit Monoclonal Antibody
产品货号: AMRe86408



This gene encodes a protein that is a member of the M49 family of metallopeptidases. This cytoplasmic protein binds a single zinc ion with its zinc-binding motif (HELLGH) and has post-proline dipeptidyl aminopeptidase activity, cleaving Xaa-Pro dipeptides from the N-termini of proteins. Increased activity of this protein is associated with endometrial and ovarian cancers. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Feb 2012]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from U-87MG, HeLa cells using AMRe86408 at 1:5000.

注意事项 (Note)

For research use only .