

产品名称: Apolipoprotein L1/APOL1 Rabbit Monoclonal Antibody
产品货号: AMRe86410



产品概述 (Summary)

产品名称 (Production Name)	Apolipoprotein L1/APOL1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Apolipoprotein L1/APOL1
别名 (Alternative Names)	APOL; APO-L; FSGS4; APOL-I
基因 ID (Gene ID)	8542
蛋白 ID (SwissProt ID)	O14791.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:44 kDa; Observed MW:44 kDa

产品名称: Apolipoprotein L1/APOL1 Rabbit Monoclonal Antibody
产品货号: AMRe86410




研究背景 (Background)

This gene encodes a secreted high density lipoprotein which binds to apolipoprotein A-I. Apolipoprotein A-I is a relatively abundant plasma protein and is the major apoprotein of HDL. It is involved in the formation of most cholesteryl esters in plasma and also promotes efflux of cholesterol from cells. This apolipoprotein L family member may play a role in lipid exchange and transport throughout the body, as well as in reverse cholesterol transport from peripheral cells to the liver. Several different transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Human plasma tissue using AMRe86410 at 1:200000.

注意事项 (Note)

For research use only .