

产品名称: CAP1 Rabbit Monoclonal Antibody
产品货号: AMRe86422



产品概述 (Summary)

产品名称 (Production Name)	CAP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CAP1
别名 (Alternative Names)	CAP; CAP1-PEN
基因 ID (Gene ID)	10487
蛋白 ID (SwissProt ID)	Q01518.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:100-1:200,ICC/IF 1:200-1:500
蛋白分子量 (Molecular Weight)	Calculated MW:52 kDa; Observed MW:52 kDa

研究背景 (Background)


产品名称: CAP1 Rabbit Monoclonal Antibody
产品货号: AMRe86422




The protein encoded by this gene is related to the *S. cerevisiae* CAP protein, which is involved in the cyclic AMP pathway. The human protein is able to interact with other molecules of the same protein, as well as with CAP2 and actin. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Aug 2016]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from K562 cells using AMRe86422 at 1:1000.



Immunofluorescence analysis of HeLa cells labelling CAP1 with AMRe86422.

注意事项 (Note)

For research use only .