

产品名称: SLC27A4/FATP4 Rabbit Monoclonal Antibody
产品货号: AMRe86426



产品概述 (Summary)

产品名称 (Production Name)	SLC27A4/FATP4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLC27A4/FATP4
别名 (Alternative Names)	IPS; FATP4; ACSVL4
基因 ID (Gene ID)	10999
蛋白 ID (SwissProt ID)	Q6P1M0.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:72 kDa; Observed MW:72 kDa

研究背景 (Background)


产品名称: SLC27A4/FATP4 Rabbit Monoclonal Antibody
产品货号: AMRe86426



This gene encodes a member of a family of fatty acid transport proteins, which are involved in translocation of long-chain fatty acids cross the plasma membrane. This protein is expressed at high levels on the apical side of mature enterocytes in the small intestine, and appears to be the principal fatty acid transporter in enterocytes. Clinical studies suggest this gene as a candidate gene for the insulin resistance syndrome. Mutations in this gene have been associated with ichthyosis prematurity syndrome. [provided by RefSeq, Apr 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from C6 cells using AMRe86426 at 1:1000.

注意事项 (Note)

For research use only .