

产品名称: Adiponectin Rabbit Monoclonal Antibody
产品货号: AMRe86430



产品概述 (Summary)

产品名称 (Production Name)	Adiponectin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Adiponectin
别名 (Alternative Names)	ACDC; ADPN; APM1; APM-1; GBP28; ACRP30; ADIPQTL1
基因 ID (Gene ID)	9370
蛋白 ID (SwissProt ID)	Q15848.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW:26 kDa; Observed MW:30 kDa

研究背景 (Background)


产品名称: Adiponectin Rabbit Monoclonal Antibody
产品货号: AMRe86430



This gene is expressed in adipose tissue exclusively. It encodes a protein with similarity to collagens X and VIII and complement factor C1q. The encoded protein circulates in the plasma and is involved with metabolic and hormonal processes. Mutations in this gene are associated with adiponectin deficiency. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Apr 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Human plasma tissue using AMRe86430 at 1:1000.

注意事项 (Note)

For research use only .