

产品名称: CD23 Rabbit Monoclonal Antibody
产品货号: AMRe86445



产品概述 (Summary)

产品名称 (Production Name)	CD23 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD23
别名 (Alternative Names)	CD23; FCE2; CD23A; IGEBF; CLEC4J; BLAST-2
基因 ID (Gene ID)	2208
蛋白 ID (SwissProt ID)	P06734.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:200,FC 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:37 kDa; Observed MW:45-50 kDa

研究背景 (Background)


产品名称: CD23 Rabbit Monoclonal Antibody
产品货号: AMRe86445




The protein encoded by this gene is a B-cell specific antigen, and a low-affinity receptor for IgE. It has essential roles in B cell growth and differentiation, and the regulation of IgE production. This protein also exists as a soluble secreted form, then functioning as a potent mitogenic growth factor. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.[provided by RefSeq, Jul 2011]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from U937 cells using AMRe86445 at 1:1000.



Western blot analysis of extracts from Raji cells using AMRe86445 at 1:1000.

注意事项 (Note)

For research use only .