

产品名称: METTL3 Rabbit Monoclonal Antibody
产品货号: AMRe86452



产品概述 (Summary)

产品名称 (Production Name)	METTL3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	METTL3
别名 (Alternative Names)	M6A; IME4; Spo8; MT-A70; hMETTL3
基因 ID (Gene ID)	56339
蛋白 ID (SwissProt ID)	Q86U44.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:500,ICC/IF 1:200-1:500,FC 1:20-1:50,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:64 kDa; Observed MW:64 kDa

研究背景 (Background)

产品名称: METTL3 Rabbit Monoclonal Antibody
产品货号: AMRe86452




This gene encodes the 70 kDa subunit of MT-A which is part of N6-adenosine-methyltransferase. This enzyme is involved in the posttranscriptional methylation of internal adenosine residues in eukaryotic mRNAs, forming N6-methyladenosine.

[provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of NIH/3T3 cells labelling METTL3 with AMRe86452.

注意事项 (Note)

For research use only .