

**产品名称: SMARCD2 Rabbit Monoclonal Antibody**  
**产品货号: AMRe86472**



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | SMARCD2 Rabbit Monoclonal Antibody     |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,IHC,FC,IP                           |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                        |

## 产品性能 (Performance)

|                     |                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                                                                                   |
| 形式 (Form)           | Liquid                                                                                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                                                     |
| 储存溶液 (Buffer)       | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt. |
| 纯化方式 (Purification) | Affinity Purification                                                                                                                                        |

## 免疫原信息 (Immunogen)

|                        |                                      |
|------------------------|--------------------------------------|
| 基因名 (Gene Name)        | SMARCD2                              |
| 别名 (Alternative Names) | SGD2; Rsc6p; BAF60B; CRACD2; PRO2451 |
| 基因 ID (Gene ID)        | 6603, 83796, 83833                   |
| 蛋白 ID (SwissProt ID)   | Q92925, Q99JR8, O54772.              |

## 产品应用 (Application)

|                          |                                                             |
|--------------------------|-------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:100-1:200,FC 1:200-1:500,IP 1:20-1:50 |
| 蛋白分子量 (Molecular Weight) | Calculated MW:59 kDa; Observed MW:59 kDa                    |

## 研究背景 (Background)


**产品名称: SMARCD2 Rabbit Monoclonal Antibody**  
**产品货号: AMRe86472**



The protein encoded by this gene is a member of the SWI/SNF family of proteins, whose members display helicase and ATPase activities and which are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The encoded protein is part of the large ATP-dependent chromatin remodeling complex SNF/SWI and has sequence similarity to the yeast Swp73 protein. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of extracts from HeLa, T24, C2C12 cells using AMRe86472 at 1:1000.

## 注意事项 (Note)

For research use only .