

产品名称: RNF14 Rabbit Monoclonal Antibody
产品货号: AMRe86475



产品概述 (Summary)

产品名称 (Production Name)	RNF14 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RNF14
别名 (Alternative Names)	ARA54; HFB30; TRIAD2; HRIHFB2038
基因 ID (Gene ID)	9604
蛋白 ID (SwissProt ID)	Q9UBS8.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,FC 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:54 kDa; Observed MW:54 kDa

研究背景 (Background)


产品名称: RNF14 Rabbit Monoclonal Antibody
产品货号: AMRe86475




The protein encoded by this gene contains a RING zinc finger, a motif known to be involved in protein-protein interactions. This protein interacts with androgen receptor (AR) and may function as a coactivator that induces AR target gene expression in prostate. A dominant negative mutant of this gene has been demonstrated to inhibit the AR-mediated growth of prostate cancer. This protein also interacts with class III ubiquitin-conjugating enzymes (E2s) and may act as a ubiquitin-ligase (E3) in the ubiquitination of certain nuclear proteins. Six alternatively spliced transcript variants encoding two distinct isoforms have been reported. [provided by RefSeq, Jan 2011]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of RNF14 in A549,HL-60,U2OS using RNF14 antibody(1:1000 diluted)



Western blot analysis of extracts from HeLa, HepG2, HCT116, T24, C2C12 cells and Rat brain tissue using AMRe86475 at 1:1000.

注意事项 (Note)

For research use only .