

**产品名称: Cpn10 Rabbit Monoclonal Antibody**  
**产品货号: AMRe86492**



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | Cpn10 Rabbit Monoclonal Antibody       |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,FC                                  |
| 种属反应性 (Reactivity)     | Human,Rat                              |

## 产品性能 (Performance)

|                     |                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                                                                                   |
| 形式 (Form)           | Liquid                                                                                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                                                     |
| 储存溶液 (Buffer)       | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt. |
| 纯化方式 (Purification) | Affinity Purification                                                                                                                                        |

## 免疫原信息 (Immunogen)

|                        |                          |
|------------------------|--------------------------|
| 基因名 (Gene Name)        | Cpn10                    |
| 别名 (Alternative Names) | EPF; CPN10; GROES; HSP10 |
| 基因 ID (Gene ID)        | 3336                     |
| 蛋白 ID (SwissProt ID)   | P61604.                  |

## 产品应用 (Application)

|                          |                                          |
|--------------------------|------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:5000,FC 1:10-1:100           |
| 蛋白分子量 (Molecular Weight) | Calculated MW:11 kDa; Observed MW:11 kDa |

## 研究背景 (Background)


**产品名称: Cpn10 Rabbit Monoclonal Antibody**  
**产品货号: AMRe86492**



This gene encodes a major heat shock protein which functions as a chaperonin. Its structure consists of a heptameric ring which binds to another heat shock protein in order to form a symmetric, functional heterodimer which enhances protein folding in an ATP-dependent manner. This gene and its co-chaperonin, HSPD1, are arranged in a head-to-head orientation on chromosome 2. Naturally occurring read-through transcription occurs between this locus and the neighboring locus MOBKL3.[provided by RefSeq, Feb 2011]

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of extracts from HepG2, HCT116 cells and Rat kidney tissue using AMRe86492 at 1:3000.

## 注意事项 (Note)

For research use only .