

产品名称: Cip4 Rabbit Monoclonal Antibody
产品货号: AMRe86493



产品概述 (Summary)

产品名称 (Production Name)	Cip4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cip4
别名 (Alternative Names)	STP; CIP4; HSTP; STOT; TRIP-10
基因 ID (Gene ID)	9322
蛋白 ID (SwissProt ID)	Q15642.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:68 kDa; Observed MW:80 kDa

研究背景 (Background)


产品名称: Cip4 Rabbit Monoclonal Antibody
产品货号: AMRe86493




Required for translocation of GLUT4 to the plasma membrane in response to insulin signaling (By similarity). Required to coordinate membrane tubulation with reorganization of the actin cytoskeleton during endocytosis. Binds to lipids such as phosphatidylinositol 4,5-bisphosphate and phosphatidylserine and promotes membrane invagination and the formation of tubules. Also promotes CDC42-induced actin polymerization by recruiting WASL/N-WASP which in turn activates the Arp2/3 complex. Actin polymerization may promote the fission of membrane tubules to form endocytic vesicles. Required for the formation of podosomes, actin-rich adhesion structures specific to monocyte-derived cells. May be required for the lysosomal retention of FASLG/FASL.

研究领域 (Research Area)

图片 (Image Data)




Western blot detection of Cip4 in 3T3 cell lysates using Cip4 antibody(1:1000 diluted).



Immunohistochemical analysis of paraffin-embedded human breast cancer using AMRe86493 antibody.

产品名称: Cip4 Rabbit Monoclonal Antibody
产品货号: AMRe86493



Western blot analysis of extracts from Rat kidney tissue using AMRe86493 at 1:3000.

注意事项 (Note)

For research use only .