

产品名称: Progesterone Receptor Rabbit Monoclonal Antibody
产品货号: AMRe86497



产品概述 (Summary)

产品名称 (Production Name)	Progesterone Receptor Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Progesterone Receptor
别名 (Alternative Names)	PR; NR3C3
基因 ID (Gene ID)	5241
蛋白 ID (SwissProt ID)	P06401.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200,ICC/IF 1:20-1:500,FC 1:200-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:99 kDa; Observed MW:90,118 kDa

产品名称: Progesterone Receptor Rabbit Monoclonal Antibody
产品货号: AMRe86497




研究背景 (Background)

This gene encodes a member of the steroid receptor superfamily. The encoded protein mediates the physiological effects of progesterone, which plays a central role in reproductive events associated with the establishment and maintenance of pregnancy. This gene uses two distinct promoters and translation start sites in the first exon to produce several transcript variants, both protein coding and non-protein coding. Two of the isoforms (A and B) are identical except for an additional 165 amino acids found in the N-terminus of isoform B and mediate their own response genes and physiologic effects with little overlap. [provided by RefSeq, Sep 2015]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from T-47D cells using AMRe86497 at 1:1000.

注意事项 (Note)

For research use only .