

产品名称: ACTL6A Rabbit Monoclonal Antibody
产品货号: AMRe86499



产品概述 (Summary)

产品名称 (Production Name)	ACTL6A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ACTL6A
别名 (Alternative Names)	Arp4; ACTL6; BAF53A; INO80K; ARPN-BETA
基因 ID (Gene ID)	86
蛋白 ID (SwissProt ID)	O96019.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:48 kDa; Observed MW:48 kDa


研究背景 (Background)

产品名称: ACTL6A Rabbit Monoclonal Antibody
产品货号: AMRe86499


This gene encodes a family member of actin-related proteins (ARPs), which share significant amino acid sequence identity to conventional actins. Both actins and ARPs have an actin fold, which is an ATP-binding cleft, as a common feature. The ARPs are involved in diverse cellular processes, including vesicular transport, spindle orientation, nuclear migration and chromatin remodeling. This gene encodes a 53 kDa subunit protein of the BAF (BRG1/brm-associated factor) complex in mammals, which is functionally related to SWI/SNF complex in *S. cerevisiae* and *Drosophila*; the latter is thought to facilitate transcriptional activation of specific genes by antagonizing chromatin-mediated transcriptional repression. Together with beta-actin, it is required for maximal ATPase activity of BRG1, and for the association of the BAF complex with chromatin/matrix. Three transcript variants that encode two different protein isoforms have been described. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of ACTL6A in HeLa, MCF-7, HepG2, U2OS using ACTL6A antibody (1:1000 diluted).



Western blot analysis of extracts from K562, Jurkat cells using AMRe86499 at 1:1000.

注意事项 (Note)

For research use only.

产品名称: ACTL6A Rabbit Monoclonal Antibody
产品货号: AMRe86499

