

产品名称: GPCR GPR49 Rabbit Monoclonal Antibody
产品货号: AMRe86510



产品概述 (Summary)

产品名称 (Production Name)	GPCR GPR49 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPCR GPR49
别名 (Alternative Names)	FEX; HG38; GPR49; GPR67; GRP49
基因 ID (Gene ID)	8549
蛋白 ID (SwissProt ID)	O75473.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:100 kDa; Observed MW:100 kDa

研究背景 (Background)


产品名称: GPCR GPR49 Rabbit Monoclonal Antibody
产品货号: AMRe86510



The protein encoded by this gene is a leucine-rich repeat-containing receptor (LGR) and member of the G protein-coupled, 7-transmembrane receptor (GPCR) superfamily. The encoded protein is a receptor for R-spondins and is involved in the canonical Wnt signaling pathway. This protein plays a role in the formation and maintenance of adult intestinal stem cells during postembryonic development. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2015]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse muscle tissue using AMRe86510 at 1:50000.

注意事项 (Note)

For research use only .