

产品名称: ITIH4 Rabbit Monoclonal Antibody
产品货号: AMRe86511



产品概述 (Summary)

产品名称 (Production Name)	ITIH4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ITIH4
别名 (Alternative Names)	H4P; IHRP; GP120; PK120; ITIHL1; PK-120; ITI-HC4
基因 ID (Gene ID)	3700
蛋白 ID (SwissProt ID)	Q14624.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:500,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:103 kDa; Observed MW:120 kDa

研究背景 (Background)


产品名称: ITIH4 Rabbit Monoclonal Antibody
产品货号: AMRe86511



The protein encoded by this gene is secreted into the blood, where it is cleaved by plasma kallikrein into two smaller forms. Expression of this gene has been detected only in liver, and it seems to be upregulated during surgical trauma. This gene is part of a cluster of similar genes on chromosome 3. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2009]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Human plasma tissue using AMRe86511 at 1:50000.

注意事项 (Note)

For research use only .