

产品名称: TOMM20 Rabbit Monoclonal Antibody
产品货号: AMRe86525



产品概述 (Summary)

产品名称 (Production Name)	TOMM20 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TOMM20
别名 (Alternative Names)	MAS20; MOM19; TOM20
基因 ID (Gene ID)	9804
蛋白 ID (SwissProt ID)	Q15388.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100,ICC/IF 1:100-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:16 kDa; Observed MW:16 kDa

研究背景 (Background)


产品名称: TOMM20 Rabbit Monoclonal Antibody
产品货号: AMRe86525




Central component of the receptor complex responsible for the recognition and translocation of cytosolically synthesized mitochondrial preproteins. Together with TOM22 functions as the transit peptide receptor at the surface of the mitochondrion outer membrane and facilitates the movement of preproteins into the TOM40 translocation pore (By similarity). [UniProt]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HepG2 cells using AMRe86525 at 1:1000.



Immunofluorescence analysis of HeLa cells labelling TOMM20 with AMRe86525.

注意事项 (Note)

For research use only .