

产品名称: BTF Rabbit Monoclonal Antibody
产品货号: AMRe86535



产品概述 (Summary)

产品名称 (Production Name)	BTF Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BTF
别名 (Alternative Names)	BTF; bK211L9.1
基因 ID (Gene ID)	9774
蛋白 ID (SwissProt ID)	Q9NYF8.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:200-1:500,ICC/IF 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW:106 kDa; Observed MW:150 kDa

研究背景 (Background)


产品名称: BTF Rabbit Monoclonal Antibody
产品货号: AMRe86535



This gene encodes a transcriptional repressor that interacts with several members of the BCL2 family of proteins. Overexpression of this protein induces apoptosis, which can be suppressed by co-expression of BCL2 proteins. The protein localizes to dot-like structures throughout the nucleus, and redistributes to a zone near the nuclear envelope in cells undergoing apoptosis. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa, 293T, RAW264.7, PC-12 cells using AMRe86535 at 1:5000.

注意事项 (Note)

For research use only .