

产品名称: USP19 Rabbit Monoclonal Antibody
产品货号: AMRe86537



产品概述 (Summary)

产品名称 (Production Name)	USP19 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	USP19
别名 (Alternative Names)	ZMYND9
基因 ID (Gene ID)	10869
蛋白 ID (SwissProt ID)	O94966.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,FC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:146 kDa; Observed MW:146 kDa

研究背景 (Background)


产品名称: USP19 Rabbit Monoclonal Antibody
产品货号: AMRe86537



Protein ubiquitination controls many intracellular processes, including cell cycle progression, transcriptional activation, and signal transduction. This dynamic process, involving ubiquitin conjugating enzymes and deubiquitinating enzymes, adds and removes ubiquitin. Deubiquitinating enzymes are cysteine proteases that specifically cleave ubiquitin from ubiquitin-conjugated protein substrates. This protein is a ubiquitin protein ligase and plays a role in muscle wasting. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2017]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from C2C12 cells using AMRe86537 at 1:1000.

注意事项 (Note)

For research use only .