

产品名称: GSDMB Rabbit Monoclonal Antibody
产品货号: AMRe86543



产品概述 (Summary)

产品名称 (Production Name)	GSDMB Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GSDMB
别名 (Alternative Names)	GSDML; PP4052; GSDMB-1; PRO2521
基因 ID (Gene ID)	55876
蛋白 ID (SwissProt ID)	Q8TAX9.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,FC 1:200-1:500,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:47 kDa; Observed MW:47 kDa

研究背景 (Background)


产品名称: GSDMB Rabbit Monoclonal Antibody
产品货号: AMRe86543



This gene encodes a member of the gasdermin-domain containing protein family. Other gasdermin-family genes are implicated in the regulation of apoptosis in epithelial cells, and are linked to cancer. Alternative splicing and the use of alternative promoters results in multiple transcript variants. Additional variants have been described, but they are candidates for nonsense-mediated mRNA decay (NMD) and are unlikely to be protein-coding. [provided by RefSeq, Nov 2016]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from A431 cells using AMRe86543 at 1:1000.

注意事项 (Note)

For research use only .