

产品名称: Cripto1 Rabbit Monoclonal Antibody
产品货号: AMRe86555



产品概述 (Summary)

产品名称 (Production Name)	Cripto1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cripto1
别名 (Alternative Names)	CR; CRGF; CRIPTO
基因 ID (Gene ID)	6997
蛋白 ID (SwissProt ID)	P13385.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000, ICC/IF 1:200-1:500, FC 1:200-1:500, IP 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:21 kDa; Observed MW:21, 18 kDa

研究背景 (Background)


产品名称: Cripto1 Rabbit Monoclonal Antibody
产品货号: AMRe86555



This gene encodes an epidermal growth factor-related protein that contains a cripto, FRL-1, and cryptic domain. The encoded protein is an extracellular, membrane-bound signaling protein that plays an essential role in embryonic development and tumor growth. Mutations in this gene are associated with forebrain defects. Pseudogenes of this gene are found on chromosomes 2, 3, 6, 8, 19 and X. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Mar 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from NCCIT cells using AMRe86555 at 1:1000.

注意事项 (Note)

For research use only .