

产品名称: PKMYT1 Rabbit Monoclonal Antibody
产品货号: AMRe86558



产品概述 (Summary)

产品名称 (Production Name)	PKMYT1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PKMYT1
别名 (Alternative Names)	MYT1; PPP1R126
基因 ID (Gene ID)	9088
蛋白 ID (SwissProt ID)	Q99640.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW:55 kDa; Observed MW:70 kDa

研究背景 (Background)


产品名称: PKMYT1 Rabbit Monoclonal Antibody
产品货号: AMRe86558



This gene encodes a member of the serine/threonine protein kinase family. The encoded protein is a membrane-associated kinase that negatively regulates the G2/M transition of the cell cycle by phosphorylating and inactivating cyclin-dependent kinase 1. The activity of the encoded protein is regulated by polo-like kinase 1. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, May 2012]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from K562 cells using AMRe86558 at 1:1000.

注意事项 (Note)

For research use only .