

产品名称: Oct-2 Rabbit Monoclonal Antibody
产品货号: AMRe86586



产品概述 (Summary)

产品名称 (Production Name)	Oct-2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	OCT2
别名 (Alternative Names)	OCT2; OTF2; Oct-2
基因 ID (Gene ID)	5452
蛋白 ID (SwissProt ID)	P09086.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000, ICC/IF 1:100-1:200, FC 1:20-1:50, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:51 kDa; Observed MW:60 kDa

研究背景 (Background)


产品名称: Oct-2 Rabbit Monoclonal Antibody
产品货号: AMRe86586



The protein encoded by this gene is a homeobox-containing transcription factor of the POU domain family. The encoded protein binds the octamer sequence 5'-ATTGTCAT-3', a common transcription factor binding site in immunoglobulin gene promoters. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Raji cells using AMRe86586 at 1:1000.

注意事项 (Note)

For research use only .