

产品名称: RBMS1 Rabbit Monoclonal Antibody
产品货号: AMRe86593



产品概述 (Summary)

产品名称 (Production Name)	RBMS1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RBMS1
别名 (Alternative Names)	YC1; MSSP-1; MSSP-2; MSSP-3; 2600014B10Rik
基因 ID (Gene ID)	5937
蛋白 ID (SwissProt ID)	P29558.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW:45 kDa; Observed MW:45 kDa

研究背景 (Background)


产品名称: RBMS1 Rabbit Monoclonal Antibody
产品货号: AMRe86593



Predicted to enable mRNA 3'-UTR binding activity; poly(A) binding activity; and poly(U) RNA binding activity. Predicted to act upstream of or within DNA replication. Predicted to be part of ribonucleoprotein complex. Predicted to be active in cytosol and nucleus. Is expressed in several structures, including alimentary system; cardiovascular system; central nervous system; genitourinary system; and immune system. Orthologous to human RBMS1 (RNA binding motif single stranded interacting protein 1). [provided by Alliance of Genome Resources, Apr 2022]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HepG2 cells using AMRe86593 at 1:1000.

注意事项 (Note)

For research use only .