

产品名称: FOXM1 Rabbit Monoclonal Antibody
产品货号: AMRe86605



产品概述 (Summary)

产品名称 (Production Name)	FOXM1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FOXM1
别名 (Alternative Names)	MPP2; HFH11; HNF-3; INS-1; MPP-2; PIG29; FKHL16; FOXM1B; HFH-11; TRIDENT; MPHOSPH2
基因 ID (Gene ID)	2305
蛋白 ID (SwissProt ID)	Q08050.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:84 kDa; Observed MW:

产品名称: FOXM1 Rabbit Monoclonal Antibody
产品货号: AMRe86605




研究背景 (Background)

The protein encoded by this gene is a transcriptional activator involved in cell proliferation. The encoded protein is phosphorylated in M phase and regulates the expression of several cell cycle genes, such as cyclin B1 and cyclin D1. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2011]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe86605 at 1:1000.

注意事项 (Note)

For research use only .