

产品名称: Integrin beta 1 Rabbit Monoclonal Antibody
产品货号: AMRe86608



产品概述 (Summary)

产品名称 (Production Name)	Integrin beta 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Integrin beta 1
别名 (Alternative Names)	CD29; FNRB; MDF2; VLAB; GPIIA; MSK12; VLA-BETA
基因 ID (Gene ID)	3688
蛋白 ID (SwissProt ID)	P05556.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:10000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:88 kDa; Observed MW:120-160 kDa

研究背景 (Background)

产品名称: Integrin beta 1 Rabbit Monoclonal Antibody
产品货号: AMRe86608




Integrins are heterodimeric proteins made up of alpha and beta subunits. At least 18 alpha and 8 beta subunits have been described in mammals. Integrin family members are membrane receptors involved in cell adhesion and recognition in a variety of processes including embryogenesis, hemostasis, tissue repair, immune response and metastatic diffusion of tumor cells. This gene encodes a beta subunit. Multiple alternatively spliced transcript variants which encode different protein isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from T24 cells using AMRe86608 at 1:1000.



Western blot analysis of extracts from Mouse small intestine, Mouse kidney tissue using AMRe86608 at 1:1000.

注意事项 (Note)

For research use only .