

产品名称: Fibrinogen alpha chain Rabbit Monoclonal Antibody
产品货号: AMRe86618



产品概述 (Summary)

产品名称 (Production Name)	Fibrinogen alpha chain Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Fibrinogen alpha chain
别名 (Alternative Names)	Fib2
基因 ID (Gene ID)	2243
蛋白 ID (SwissProt ID)	P02671.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:500-1:2000,ICC/IF 1:20-1:50,FC 1:20-1:50,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:95 kDa; Observed MW:110 kDa

产品名称: Fibrinogen alpha chain Rabbit Monoclonal Antibody
产品货号: AMRe86618




研究背景 (Background)

This gene encodes the alpha subunit of the coagulation factor fibrinogen, which is a component of the blood clot. Following vascular injury, the encoded preproprotein is proteolytically processed by thrombin during the conversion of fibrinogen to fibrin. Mutations in this gene lead to several disorders, including dysfibrinogenemia, hypofibrinogenemia, afibrinogenemia and renal amyloidosis. Alternative splicing results in multiple transcript variants, at least one of which encodes an isoform that undergoes proteolytic processing. [provided by RefSeq, Jan 2016]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HepG2 cells using AMRe86618 at 1:1000.

注意事项 (Note)

For research use only .