

产品名称: Integrin alpha 2 Rabbit Monoclonal Antibody
产品货号: AMRe86622



产品概述 (Summary)

产品名称 (Production Name)	Integrin alpha 2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Integrin alpha 2
别名 (Alternative Names)	BR; GPIa; CD49B; HPA-5; VLA-2; VLAA2
基因 ID (Gene ID)	3673
蛋白 ID (SwissProt ID)	P17301.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:129 kDa; Observed MW:150 kDa

研究背景 (Background)


产品名称: Integrin alpha 2 Rabbit Monoclonal Antibody
产品货号: AMRe86622



This gene encodes the alpha subunit of a transmembrane receptor for collagens and related proteins. The encoded protein forms a heterodimer with a beta subunit and mediates the adhesion of platelets and other cell types to the extracellular matrix. Loss of the encoded protein is associated with bleeding disorder platelet-type 9. Antibodies against this protein are found in several immune disorders, including neonatal alloimmune thrombocytopenia. This gene is located adjacent to a related alpha subunit gene. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from A549 cells using AMRe86622 at 1:1000.

注意事项 (Note)

For research use only .