

产品名称: eNOS Rabbit Monoclonal Antibody
产品货号: AMRe86625



产品概述 (Summary)

产品名称 (Production Name)	eNOS Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	eNOS
别名 (Alternative Names)	eNOS; ECNOS
基因 ID (Gene ID)	4846
蛋白 ID (SwissProt ID)	P29474.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:133 kDa; Observed MW:140 kDa

研究背景 (Background)


产品名称: eNOS Rabbit Monoclonal Antibody
产品货号: AMRe86625




Nitric oxide is a reactive free radical which acts as a biologic mediator in several processes, including neurotransmission and antimicrobial and antitumoral activities. Nitric oxide is synthesized from L-arginine by nitric oxide synthases. Variations in this gene are associated with susceptibility to coronary spasm. Alternative splicing and the use of alternative promoters results in multiple transcript variants. [provided by RefSeq, Oct 2016]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Rat heart tissue using AMRe86625 at 1:1000.



Western blot analysis of extracts from Mouse spleen tissue using AMRe86625 at 1:1000.

注意事项 (Note)

For research use only .